Swemed™ Transducer Covers

26226.01

EN: Instructions for use – Swemed™ Transducer Covers

Product Description

Transducer Cover is a sterile, latex-free cover packed with two clips. It is intended for covering the ultrasound transducer during follicle aspiration procedures.

Instructions for use

Check that the package is intact and the "use by" date has not expired.

Carefully open the bag containing the transducer cover.

Place the transducer cover over the ultrasound transducer used for vaginal ultrasound procedures.

Use the cable clips to keep the cover in place.

If Transducer Cover is to be used, apply the gel over the ultrasound transducer head and insert the transducer in a cover.

Precautions

Discard the product if the packaging or product is damaged.

The product may only be used by, or under the guidance of, a person qualified in gynaecological and follicle aspiration procedures.

Caution: Federal (US) law restricts this device to sale by, or on the order of, a physician.

Specifications

STERILE	Sterilized using irradiation.
	Bacterial Endotoxins (LAL assay) ≤ 1.0 EU/device.
	Do not re-use, discard after procedure.
	Use by – see label.
	Caution: Consult accompanying documents.
REF	Catalog number
LOT	Batch code.

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Vitrolife Sweden AB, Box 9080 (Gustaf Werners gata 2), SE-400 92 Göteborg, Sweden.
Vitrolife Inc., 3601 South Inca Street, Englewood CO 80110, USA.
Vitrolife Pty. Ltd., 55 Flemington Rd, North Melbourne VIC 3051, Australia.
Vitrolife Sweden AB